

HPLC API调试

API调试记录 by 孙工

1. 获取所有项目

HTTP

My Collection > Post data

Save

Share

POST

10.10.10.1:8001/project/findAllProjectList

Send

Docs

Params

Authorization

Headers 9

Body

Scripts

Tests

Settings

Cookies

Headers

8 hidden

Key	Value	Descri...	Bulk Edit	Presets	...
☒ Content-Type	application/json;charset=UTF-8				
Key	Value	Description			

Body

200 OK • 2.21 s • 301 B • Save Response

JSON

Preview

Visualize

1

2

3

4

5

6

7

8

9

10

11

12

13

14

{

"code": "0",

"errorCode": null,

"msg": "success",

"data": [

{

"projectId": 1,

"projectName": "default-project"

},

{

"projectId": 2,

"projectName": "Test"

}

]

}

2. 获取项目中的序列方法

HTTP My Collection > Post data Save Share

POST 10.10.10.1:8001/project/findSequenceMethodList Send

Docs Params Authorization Headers 9 Body Scripts Tests Settings Cookies

raw JSON Schema Beautify

```
1 {
2   |   "projectId": 2
3   }
```

Body 200 OK • 21 ms • 244 B • Save Response

{ JSON Preview Visualize

```
1 {
2   |   "code": "0",
3   |   "errCode": null,
4   |   "msg": "success",
5   |   "data": [
6   |     {
7   |       |   "id": 5,
8   |       |   "methodName": "seq"
9   |     }
10  |   ]
11 }
```

3. 获取项目中的报告方法

POST 10.10.10.1:8001/project/findReportMethodList Send

Docs Params Authorization Headers 9 Body Scripts Tests Settings Cookies

none form-data x-www-form-urlencoded raw binary GraphQL JSON Schema Beautify

```
1 {
2   "projectId": 2
3 }
```

Body Cookies Headers 5 Test Results 1/1 200 OK • 121 ms • 783 B Save Response

{ JSON Preview Visualize

```
1 {
2   "code": "0",
3   "errCode": null,
4   "msg": "success",
5   "data": [
6     {
7       "id": 27,
8       "methodName": "结果综合报告-1"
9     },
10    {
11      "id": 28,
12      "methodName": "3D匹配图报告"
13    },
14    {
15      "id": 29,
16      "methodName": "3D等高图报告"
17    },
18    {
19      "id": 17,
20      "methodName": "单结果报告-2"
21    },
22    {
23      "id": 18,
24      "methodName": "样品综合报告"
25    },
26  ]
27 }
```

4. 获取项目中的处理方法

POST 10.10.10.1:8001/project/findProcessMethodList Send

Docs Params Authorization Headers 9 Body Scripts Tests Settings Cookies

none form-data x-www-form-urlencoded raw binary GraphQL JSON Schema Beautify

```
1 {
2   "projectId": 2
3 }
```

Body Cookies Headers 5 Test Results 1/1 200 OK • 13 ms • 250 B Save Response

{ JSON Preview Visualize

```
1 {
2   "code": "0",
3   "errCode": null,
4   "msg": "success",
5   "data": [
6     {
7       "id": 32,
8       "methodName": "process1"
9     }
10  ]
11 }
```

5. 获取仪器系统

需要设备在线才能获取到

POST 10.10.10.1:8001/project/findInstrumentList Send

Docs Params Authorization Headers 8 Body Scripts Tests Settings Cookies

none form-data x-www-form-urlencoded raw binary GraphQL JSON Schema Beautify

```
1
```

Body Cookies Headers 5 Test Results 1/1 200 OK • 1.18 s • 266 B Save Response

{ JSON Preview Visualize

```
1 {
2   "code": "0",
3   "errCode": null,
4   "msg": "success",
5   "data": [
6     {
7       "id": 1,
8       "instrumentName": "HPLC",
9       "onlineStatus": 1
10    }
11  ]
12 }
```

6. 新建序列方法

注：每次新建序列方法都要生成一个新的流程Id：procInstId。后续运行序列、控制序列、查询序列状态都需要用同样的流程Id进行关联，否则后续操作会失败

POST10.10.10.1:8001/project/addSequenceMethodSend

DocsParamsAuthorizationHeaders9BodyScriptsTestsSettingsCookies

noneform-datax-www-form-urlencodedrawbinaryGraphQLJSONSchemaBeautify

```
1 {
2   "procInstId": "3",
3   "projectId": 2,
4   "sequenceMethodId": 5,
5   "name": "seq2",
6   "seqTrayNum": 1,
7   "instIdArr": [1],
8   "sampleArr": [
9     {
10      "trayNum": 0,
11      "xPosition": "A",
12      "yPosition": 1,
13      "sampleName": "s1"
14    }
15  ]
16 }
```

BodyCookiesHeaders5Test Results1/1200 OK95 ms219 BSave Response

{ }JSONPreviewVisualize

```
1 {
2   "code": "0",
3   "errCode": null,
4   "msg": "success",
5   "data": null
6 }
```

7. 运行序列

流程id必须与创建序列方法一致

POST 10.10.10.1:8001/project/runSequence Send

Docs Params Authorization Headers 9 Body Scripts Tests Settings Cookies

none form-data x-www-form-urlencoded raw binary GraphQL JSON Schema Beautify

```
1 {
2   "procInstId": "3",
3   "proMethodId": 32,
4   "reportMethodId": 27,
5   "shutDown": 0,
6   "exportPath": "D:\\\"
7 }
8
```

Body Cookies Headers 5 Test Results 1/1 200 OK • 1.47 s • 219 B • Save Response

{ JSON Preview Visualize

```
1 {
2   "code": "0",
3   "errCode": null,
4   "msg": "success",
5   "data": null
6 }
```

8. 序列控制

POST 10.10.10.1:8001/project/controlSequence Send

Docs Params Authorization Headers 9 Body Scripts Tests Settings Cookies

none form-data x-www-form-urlencoded raw binary GraphQL JSON Schema Beautify

```
1 {
2   "action":1
3 }
4
```

Body Cookies Headers 5 Test Results 1/1 200 OK • 35 ms • 219 B • Save Response

{ } JSON Preview Visualize

```
1 {
2   "code": "0",
3   "errCode": null,
4   "msg": "success",
5   "data": null
6 }
```

9. 获取序列状态

流程id必须与创建序列方法一致

POST 10.10.10.1:8001/project/getSequenceStatus Send

Docs Params Authorization Headers 9 Body Scripts Tests Settings Cookies

none form-data x-www-form-urlencoded raw binary GraphQL JSON Schema Beautify

```
1 {
2   |   "procInstId":3
3   | }
4 }
```

Body Cookies Headers 5 Test Results 1/1 200 OK • 15 ms • 242 B | Save Response

{ } JSON Preview Visualize

```
1 {
2   |   "code": "0",
3   |   "errCode": null,
4   |   "msg": "success",
5   |   "data": {
6   |     |   "seqTrayNum": 1,
7   |     |   "status": 1
8   |   | }
9   | }
```

10. 序列运行完成后，D盘能够生成报告